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Introduction

Thank you for purchasing a Caremaxx doppler. This doppler is a quality product from Caremaxx.

Carefully read this user manual, before you start using the device. Keep the user manual for future reference.

Safety precautions

- Only ever use the product as described in the user manual.
- In case of improper use, the warranty will become invalid.
- The product is not intended for commercial use.
- Repairs may only be performed by authorized service workshops. Do not open the device yourself.
- Keep the device out of reach from babies, children or intellectually incompetent persons.
- Make sure your hands are dry when you pull the plug out of the socket.
- The device and its accessories are not suitable for the dishwasher.

The package contains following components:

- 1x Doppler unit
- 1x headphone
- 1x Audio cable
- 1x 9V battery
- 1x manual

The packaging material can be reused or recycled. Make sure that the used packaging material ends up in the right waste processing system. If, while removing the packaging, you notice transport damage, please contact your supplier immediately.

The Caremaxx doppler is a single-hand device for detection of the fetal heart. The fetal heart may be detected from early pregnancy, and more easily from pregnancies greater than 12 weeks. The doppler is suitable for use in the antenatal.

Every doppler is supplied with 9V battery. Dry-cell battery must NOT be placed in a re-charging device and must NOT be disposed of by burning. The doppler operates at a nominal frequency of 3.3MHz, using Doppler ultrasound. There has an 'on/off/volume key(1)', it's easy to operate. And working indicator light (2) shows working condition. The doppler is lightweight and designed to fit comfortably in hand.

Operation Manual

- Install the battery correctly, turn on the power knob, make sure the light on.
- Before listening, make the volume to medium to protect your ears from unpleasant voice.
- Put the earphones on, apply the gel and slowly move the doppler until you hear your baby's heartbeat. Place the transducer on the examination site and move it around slowly until a good signal is obtained. This will sound like 'galloping horses'. It may be necessary to angle the transducer slightly to obtain the optimum signal.
- You can adjust the volume knob to catch the optimum signal.
- Fetal heart sounds is different from yours, fetal heart rate is 120-160bpm, it sounds like drumming. Please pay your patience to listen, you will get it.
- It's possible to hear sounds from 16 weeks into pregnancy.

The best period of time to hear Fetal Heart Rate movement are:

1. Within 30 minutes after the pregnant woman gets up in the morning.
2. Within 60 minutes after the pregnant woman has lunch.
3. 30 minutes before pregnant woman goes to bed at night.

Use of audio line

Insert one jack of the audio line into the main unit's headset socket (anyone of the two), put another plug into the recorder, then the fetal heart beat can be record by recorder.

Move the transducer to pregnant woman's left chest to record her heart

beat and play the mother's heart beat to the baby after the baby is born, this can calm the baby.

To switch off turn 'on/off/volume key(1)' .

Battery fitting or replacement

To fit or replace the battery on the doppler:

Remove the rear battery compartment cover(5) carefully . Remove the old battery. Fit the new battery (IEC6F22 9V alkaline). Then off the battery cover.

Use only the specified battery type and insert the battery according to polarity instruction.

Always dispose of empty batteries in accordance with regulations. Do not dispose together with household garbage

New alkaline battery will give more than 500 one minute examinations. Battery should be removed if the unit is not in regular use.

Preventive maintenance

After use, the doppler should be dried thoroughly with a soft cloth. If the doppler is soiled after use, follow the cleaning, and disinfection procedure detailed in the preventive maintenance section.

General

The equipment is designed to require a minimum amount of maintenance. To obtain the best performance and maintain safety, the following checks should be carried out quarterly or annually, depending on usage.

Check the doppler for damage or cracks which may allow the ingress of liquids or gel.

Cleaning and disinfecting

Cleaning: Wipe the instrument case with a cloth dampened in soap or a detergent solution and wipe dry with a clean cloth.

Disinfection: If soiled, clean as above, then wipe the instrument case with an alcohol-impregnated(70% ethanol or isopropyl).

Notes:

The doppler prenatal heart listener brings you a clear sounds of your baby's heartbeat , but in order to protect your ears from unpleasant voice, we suggest you power on the device first before listening and make the volume in middle place.

Your baby and your body will make different sounds, at the time to listen, you can adjust the volume for the best effect.

You will be suggested to pay some patience to find the best position of fetal heartbeat. Heartbeat sounds level has relations with fetal heart position. Sometimes, you may not hear because of the baby's movement.

In order to obtain a good sound effect, you're suggested to use Gel. If you do not use Gel, you also can put food oil or water. The doppler probe should be touched your skin. Mobile phone, radio signals, etc will effect the Doppler's use. Please keep them away you. We suggest you not to detect fetal heartbeat more than 3 times a day, each time less than 3 minutes.

Trouble Shooting

If the doppler do not work, please check follows:

1. Volume knob whether or not put the lowest
2. Bad or inadequate battery power
3. Earphone cable whether or not good contact with the doppler's Earphone jack.
4. Whether or not put the Gel
5. If you see the light on but the doppler was not working, that's lower battery caused please change a new battery.

Notes: it's a normal electricity sounds when you turn on the doppler, not a bad or headphone's problems.

Disclaimer

The doppler is a non diagnostic device and does not replace the role of the doctor. You ought to consult the doctor about any pregnancy question. However, we still try our best to offer you the fetal heart rate for your reference. Do not as according as clinic diagnose but reference only during gestation. At the same time please operation it as instruction manual. Our industry and dealer will carefree if any of inadequately use.

General

Diagnostic ultrasound has been in use for over 25 years with no confirmed adverse effects on patients or instrument operators at the intensities typical of present diagnostic instruments.

Although the total absence of adverse effects to human subjects after such extensive use at diagnostic power levels is gratifying, available data are not conclusive and the possibility that biological effects may be identified in the future remains.

It is therefore deemed desirable by medical and other scientific authorities in this field that exposure to ultrasound should be limited to the duration and intensity appropriate for the clinical objective. Because fetal tissue could be more sensitive to biological effects by reason of pregnant subjects be kept to a minimum.

At present, there is a clear consensus that the benefits to patients of prudent use of diagnostic ultrasound outweigh the risks, if any, that may be present.

The Caremaxx Doppler is a portable battery operated detector designed for the detection of fetal life and confirmation of continued life during pregnancy.

Ultrasound safety considerations and data

Minimizing patient exposure

Acoustic output of the doppler is internally controlled and cannot be changed by the operator in the course of the examination. The duration of use is,

however, fully under the control of the operator. Mastery of the techniques described in the operating instructions will facilitate limit the maximum amount of diagnostic information with the minimum of exposure.

Acoustic output data

The acoustic output of the doppler transducer has been measured in water using a calibrated hydrophone at Authorized Laboratory. Normalized values, which estimate the maximum 'in-situ' dosage to tissue at the point of highest intensity in the beam path have been calculated.

Waste management

This instrument may not disposed of with the household waste. Consumers are obliged to send all electric or electronic instruments, whether or not they contain any hazardous materials, to a depository for environmental pollutants in their home town or to their specialized dealership to ensure their environmentally sound disposal. If you have any questions regarding waste management, please contact your municipality or your specialized vendor.

Technical specifications

Name:	Caremaxx doppler
Item number:	31021
Dimensions (D x W x H):	17 x 15 x 6.8 cm
Weight:	approx. 1 kg
Battery:	9V
EAN code:	8 717964 31021 3
Temperature:	-10°C tot 60°C
Humidity:	0% tot 95%

This CE mark on a product denotes conformity with the European Council Directive 93/42/EEC(MDD) concerning medical devices.

Warranty and conditions for repair

In cases where you want to use your warranty, you must always contact your dealer or our customer service. If you need to send the device to our customer service, you must include the complaint and a copy of your receipt.

The following conditions apply to the warranty:

1. For Caremaxx products, a warranty period applies of three years after the purchase date. This can be proven by means of the receipt or invoice.
2. All complaints resulting from defective material or manufacturing faults, will be repaired for free within the warranty period.
3. Cases where warranty applies, do not automatically lead to extension of the warranty period, neither for the device itself, nor for the parts that have been replaced.
4. Excluded from the warranty are:
 - a. All damages caused by incompetent handling, e.g. not correctly following the user manual.
 - b. Damages caused by repairs performed by the buyer or another unauthorized person.
 - c. Transport damages, caused during transport from the seller to the user or when sending the device to the customer service.
 - d. Accessories that are subject to wear.
5. The manufacturer does not accept any liability for direct or indirect consequential damage caused by the device. This also applies, if the damage to the device is recognized under this warranty.

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