Caremaxx®

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E-mail : info@caremaxx.nlSite Internet : www.caremaxx.nl

1 Caremaxx baby scale

1.1 Thank you for purchasing a Caremaxx baby scale.

Read this User Manual carefully before using the baby scale.

1.2 Safety measures

- This baby scale should strictly be used in the manner described in this User Manual.
- The guarantee will lapse in case of improper use.
- The product is not intended for commercial use.
- Repairs may only be carried out at authorised service workshops.
- The appliance should not be exposed to extremely high or low temperatures or direct sunlight.
- Place the baby scale on a flat and stable base surface (for example a commode) before use.
- Keep the weighing scale clean.
- Do not load the weighing scale with more than 20 kg.

1.3 Delivery and packaging

We want our high-quality products to reach you in the condition in which we sent them to you. That is why your Caremaxx baby scale has been packed in a stable, protected packaging made of environment-friendly materials. If you do not wish to use the box to store the appliance, you can dispose of it along with your domestic waste paper. Contact your supplier immediately if you discover transport damage when unpacking.

The packaging contains the following components:

- 1 Caremaxx baby scale, model 50300
- 4 Batteries (AA, 1.5V)
- 1 User Manual

2 Starting to use your baby scale

2.1 Inserting and removing batteries.

Turn the weighing scale over.

Open the cover of the battery compartment and insert 4 AA-batteries with the plus (+) and the minus (-) in the manner indicated in the battery compartment. Replace the cover and click it shut.

Replace the batteries as soon as the display reads "L o w".

Remove the batteries if you are not going to use the baby scale for some time, in order to prevent battery leakage into the appliance.

2.2 Functions

ON/TARE

Press the ON/TARE button to switch on the weighing scale.

The scale will now read 0.000 and is ready to take measurements.

Place the baby in the centre of the weighing scale and read the weight on the display.

You may place the baby on a thin mattress or towel without including the weight of such mattress or towel. Switch on the weighing scale. Now place the towel on the weighing scale and then press the ON/TARE button.

The weighing scale will now assume the initial position of 0.000. Now place the baby on the weighing scale and read the weight on the display.

Switch off the weighing scale by keeping the ON/TARE button pressed for 3 seconds. If you do not switch off the weighing scale, it will automatically switch off after 5 minutes.

MODE

Using the MODE button, you can change the display from kgs to lbs and vice versa. You can see which scale is active in the right part of the display.

MEM

The weight displayed can be saved in the memory by pressing the MEM button. When the weighing scale shows a reading of "0.000", you can recall the last-saved weight by pressing the MEM button.

The weight remains in the memory until you save a new weight or until the battery is discharged.

3 Other

3.1 Cleaning and maintenance

Switch off the baby scale before cleaning it. Clean it with a moist cloth, but avoid the entry of water into the appliance. Do not use any aggressive cleaning products.

Store the weighing scale in a dry place, outside direct sunlight. Remove the batteries from the appliance if you do not

intend to use the scale for some time. This will prevent the battery from leaking into the appliance.

3.2 Waste management

This instrument is not to be disposed of with the household waste. Consumers are obliged to send all electric or electronic instruments, whether or not they contain any hazardous materials, to a depository for environmental pollutants in their home town or to their specialized dealership to ensure their environmentally sound disposal. Do not forget to remove the battery from the instrument before throwing it away. Used batteries do not belong into the household waste, but should be taken to the pertinent waste management facility or a special collecting point at the supermarket or the electronics store. If you have any questions regarding waste management, please contact your municipality or your specialized vendor.

3.3 Warranty and repair conditions

If you wish to submit a warranty claim, please contact your supplier or our customer service. If the instrument is to be sent to our customer service, make sure to enclose a description of the complaint as well as a copy of your receipt of purchase.

Warranty is granted subject to the following conditions:

1. The Caremaxx products are subject to a warranty period of three years, starting with the purchasing date as evidenced by the receipt or invoice.

2. All complaints due to material or manufacturing defects will be remedied free of charge within the warranty period.
3. A warranty case does not automatically imply an extension of the warranty period; neither for the instrument itself nor for the replaced parts.
4. Excluded from the warranty are:
 - a. Any damages caused by inappropriate use; for instance due to noncompliance with the user manual.
 - b. Damages caused by repair jobs performed by the purchaser or by any other unauthorized person.
 - c. Transport damages in transit accrued enroute from the seller to the consumer or from the consumer to the customer service.
 - d. Accessories that are subject to wear.
5. The producer accepts no liability for direct or indirect consequential damages caused by the instrument; not even if the damages to the instrument constitute an acknowledged warranty case.

3.4 Technical Data

Product	Caremaxx baby scale
Model	50300
Maximum weighing capacity	20kg
Weight indication	per 5g
Batteries	4x AA (1.5V)
Dimensions:	60 x 46 x 18 cm
EAN-Code	87 17964 50300 4

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