# Caremaxx®

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Infrarood massage apparaat

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Site Internet : [www.caremaxx.nl](http://www.caremaxx.nl)

# 1 Caremaxx infra-red massage device

## 1.1 Thank you for purchasing a Caremaxx infra-red massage device.

Please carefully read these user instructions before using the infra-red massage device.

## 1.2 Safety measures

- Only use this massage device for the purposes described in this user manual.
- The guarantee becomes void in the event of improper use.
- The massage device is only suitable for dry use.
- Do not use this massage device in moist places.
- Do not use this device with wet hands.
- Use this massage device on the body, and not the hands.
- This device is not intended for commercial use.
- Do not leave this device unsupervised within the reach of children.
- Repairs may only be carried out at authorised service points.
- Do not expose the device to extremely high or low temperatures or direct sunlight.
- Do not dip any part of the massage device in water or any other liquid.
- Do not use this device under blankets or cushions.
- Do not use the device for more than 20 minutes. Allow the device to cool down for 30 minutes after use.

## Delivery and packaging

We want our quality products to reach you in the same condition that we sent them. That

is why your Caremaxx massage device comes in strong protective packaging which is made from environment-friendly materials. If you do not want to use the box to store the device, then you can simply dispose of it like normal paper. If you discover transport damage while unpacking, then please contact your supplier.

The package contains following components:  
1x Caremaxx infra-red massage device, model 30100  
1x User manual

## 2 Using your massage device

### 2.1 Massage operating positions

Ensure that the On/Off switch ③ is in the Off position before you plug the device into the mains. Then place the switch ③ into the desired position.

The massage device has 4 different operating positions:

- 1) The device is off.
- 2) Massage only
- 3) Only infra-red heat
- 4) Massage and infra-red heat at the same time

Turn off the device and remove the plug from the mains after use.

### 2.2 Usage

You can massage yourself or your partner. Massage as long as it feels comfortable and relaxing. In general, the optimum massage time lies between 10 to 15 minutes. The treatment time for infra-red heat is also 10 to 15 minutes. If you massage too long, this could lead to tense muscles instead of relaxed muscles.

Do not use the device too long on the same place and make circular movements. Hold the massage head softly against the body and avoid using too much pressure.

The device may heat up if the infra-red feature is being used. Thus do not use the device on synthetic materials, flammable or heat-sensitive materials. Allow the device to cool down before storing.

The massage device operates on the basis of vibration massage. This involves the massage head being vibrated at high frequency, with this movement being in-turn transferred to the body. If desired, infra-red heat can also be added. It is also possible to use infra-red heat separately. This special heat-application method helps to improve the massage effect.

The unique combination of intense massage and heat-radiation leads to great relaxation.

### **3 Warnings**

This massage device is not a medical device. A massage treatment at home cannot replace medical therapy.

Do not use the device on:

- The head or the face;
- Bones or joints;
- The front of the neck;
- Genitals;
- Kidneys.

Do not use the device in the following cases:

- During pregnancy;
- If you have a pacemaker, artificial joints or electronic implants;

- If you suffer from one or more of the following conditions or illnesses:
  - Thrombosis, blood circulation problems
  - Varicose veins
  - Open wounds, bruises,
  - Irritable skin or burn wounds
  - Infections

Stop the massage treatment if it feels uncomfortable or causes pain, and contact your doctor.

## 4 General

### 4.1 Cleaning and maintenance

Switch off the device and remove the plug from the mains before you clean the device. Allow the device to cool down before cleaning. Clean it using a dry or slightly moist cloth. Do not use aggressive substances, such as scouring and cleaning products, because they may damage the surface. Never dip the device in water and avoid moisture from entering it. The massage device must be completely dry before it can be used again. Store the massage device in a dry place, away from direct sunlight. Neatly roll up the power cord in order to prevent damage.

### 4.2 Waste management

This instrument is not to be disposed of with the household waste. Consumers are obliged to send all electric or electronic instruments, whether or not they contain any hazardous materials, to a depository for environmental pollutants in their home town or to their specialized dealership to ensure their environmentally sound disposal. If you have any questions regarding waste

management, please contact your municipality or your specialized vendor.

### **4.3 Warranty and repair conditions**

If you wish to submit a warranty claim, please contact your supplier or our customer service. If the instrument is to be sent to our customer service, make sure to enclose a description of the complaint as well as a copy of your receipt of purchase.

Warranty is granted subject to the following conditions:

1. The Caremaxx products are subject to a warranty period of three years, starting with the purchasing date as evidenced by the receipt or invoice.
2. All complaints due to material or manufacturing defects will be remedied free of charge within the warranty period.
3. A warranty case does not automatically imply an extension of the warranty period; neither for the instrument itself nor for the replaced parts.
4. Excluded from the warranty are:
  - a. Any damages caused by inappropriate use; for instance due to noncompliance with the user manual.
  - b. Damages caused by repair jobs performed by the purchaser or by any other unauthorized person.
  - c. Transport damages in transit accrued enroute from the seller to the consumer or from the consumer to the customer service.
  - d. Accessories that are subject to wear.
5. The producer accepts no liability for direct or indirect consequential damages caused by the instrument; not even if the damages to the instrument constitute an acknowledged warranty case.

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